

Supplemental Materials




Figure S1: 3D printed sensor block.

The 3D printed sensor block used to house the sensors of the RAMOS device.




Figure S2: Raw oxygen sensor value during measurement phases.

I: Equilibration phase, II: Measurement phase, III: Flush phase. Progress of the oxygen partial pressure during a measurement phase. The O₂ content decreases up to 10 % during the measurement phase, resulting in a dynamically changing O₂ gradient.




Figure S3: Raw oxygen sensor value of RAMOS experiment.

Progress of the oxygen partial pressure during the cultivation of Figure 3 with decreasing oxygen content for the first 12 to 14 hours.




Figure S4: Effects of CO₂, H₂ and humidity on CH₄ sensor reading.

CH₄, CO₂, and H₂ sensor signals for different gas composition and humidity. 1. dry nitrogen gas; 2. humidified nitrogen gas; 3. CO₂ ramp from 0 – 20 vol%; 4. H₂ ramp from 0 – 80 vol%.




Figure S5: Single line data of Fig. 4.

Single line data of Fig. 4. Same color/icons symbolize same process parameters. Dashed lines indicate duplicates.




Figure S6: CQ of *M. capsulatus* cultivation.

CQ of the cultivation shifts from around 0.5 during the exponential growth phase towards 0.7 during the transition phase into the O₂ limitation. During the O₂ limitation, a CQ of around 0.9 can be measured.




Figure S7: Comparison of transfer rates of *M. capsulatus* under various KNO₃ concentrations.

(a) OTR and (b) MTR courses of *M. capsulatus* under various KNO₃ concentrations. The cultivations were conducted in duplicates with the min/max values illustrated as shadows. Due to sensor errors, the 1.0 g/L and 0.15 g/L KNO₃ culture MTR data are only shown as singular courses. Before inoculation, washing the cells removed residual KNO₃ from the preculture. Cultivation conditions: NMS medium, c_N: 0 g/L – 1 g/L KNO₃, c_{buffer} = 30 mM MOPS, T = 37 °C, n = 350 rpm, d₀ = 50 mm, V_L = 10 mL, initial OD = 0.1, initial pH = 6.8, N = 2, Only every 20th data point over time is indicated by the corresponding symbol;

c_{in-gas} = 5.9 vol% CH₄, 3.4 vol% O₂, 1.3 vol% CO₂, 89.4 vol% N₂.




Figure S8: Integrals of transfer rates *M. capsulatus* under various KNO₃ concentrations.

Integrals of (a) OTR and (b) MTR courses of *M. capsulatus* under various KNO₃ concentrations. The time of the shifts for the 0.15 g/L KNO₃ and 0.3 g/L KNO₃ cultures are indicated by the vertical dashed lines. Cultivation conditions: NMS medium, c_N: 0 g/L – 1 g/L KNO₃, c_{buffer} = 30 mM MOPS, T = 37 °C, n = 350 rpm, d₀ = 50 mm, V_L = 10 mL, initial OD = 0.1, initial pH = 6.8, N = 2, Only every 20th data point over time is indicated by the corresponding symbol;

c_{in-gas} = 5.9 vol% CH₄, 3.4 vol% O₂, 1.3 vol% CO₂, 89.4 vol% N₂.




Figure S9: OQ of *M. capsulatus* cultivation under various KNO₃ concentrations.

OQs of the cultivation first 72 hours of the cultivation are shown. Both 1 g/L and 0.3 g/L KNO₃ start with an OQ of 1.0 and settle around 1.5 before the oxygen limitation. Cultures with 1 g/L KNO₃ reaches the OQ of 1.9 after 24 hours, while 0.3 g/L KNO₃ reaches the same OQ after 15 hours shortly after the nitrogen shift. The cultures with 0.15 g/L KNO₃ start around 1.45 and reach 1.9 after 13.5 hours. If no KNO₃ is present in the medium, the OQ starts at 1.9 from the beginning once the calculation threshold of OTR > 2 mmol/L/h is reached. Ratios with OTR values above 2.0 are shown. Only every 20th data point over time is indicated by the corresponding symbol.




Figure S10: CQ of *M. capsulatus* cultivation under various KNO₃ concentrations.

CQs of the cultivation first 72 hours of the cultivation are shown. Both 1 g/L and 0.3 g/L KNO₃ start with a CQ of 0.5 to 0.6 and settle around 0.7 to 0.8 after the oxygen limitation. Cultures with 1 g/L KNO₃ reaches the CQ of 0.9 after 24 hours, while 0.3 g/L KNO₃ reaches the same CQ after 15 hours shortly after the nitrogen shift. The cultures with 0.15 g/L KNO₃ start around 0.6 and reach 0.8 after 13.5 hours. If no KNO₃ is present in the medium, the CQ starts at 0.8 from the beginning once the calculation threshold of CTR > 1 mmol/L/h is reached. Ratios with CTR values above 1.0 are shown. Only every 20th data point over time is indicated by the corresponding symbol.